



THE HARD-CODED KEY TO MY HEART - HACKING A PACEMAKER PROGRAMMER


ANDERS B. WILHELMSEN @anderbw

EIVIND S. KRISTIANSEN @Skjelmo

MARIE MOE @MarieGMoe


THE STORY OF MY HEART

WHY WE STARTED A HACKING PROJECT



- Proprietary software
- No transparency
- Legacy technology
- Added connectivity = added attack surface

EMBODIED AND CONNECTED



BIOTRONIK

BITFLIP IN THE AIR GAVE US ACCESS TO A CRASH FILE FROM MY PACEMAKER

BOEING[®]
PHANTOM WORKS

Neutron Environment in the Atmosphere

Boeing Radiation Effects Lab

Neutrons, created by cosmic ray interactions with the O₂ and N₂ in the air, peak at ~60,000 ft. At 30,000 ft the neutrons are about 1/3 the peak flux, and on the ground, ~1/400 of the peak flux. The peak flux is ~4 neutron/cm²sec. Other particles such as secondary protons and pions are also created, but for SEU the neutrons are the most important.


RADIATION ENVIRONMENT IN THE ATMOSPHERE

Legend:

- Y Primary Cosmic Rays
- A Secondary Cosmic Rays (muons, pions)
- Neutrons

5








DATA ERROR IN PACEMAKER

ENCRYPTED DATA FROM MY HEART

0 0070:	30 78 30 30 30 30 30 30 46	0D 0A 38 31 20 36 64 20	0x00000F ..81 60
0 0080:	30 30 20 30 30 20 30 30 30	20 30 30 20 30 33 20 32	00 00 00 00 03
0 0090:	34 20 30 30 20 30 30 20	38 30 20 30 30 20 30 30	4 00 00 80 00 00
0 00A0:	20 31 38 20 30 30 20 30	30 20 2A 20 30 78 30 30	18 00 00 0 * 0x0
0 00B0:	30 30 30 30 20 20 30 30	30 30 20 30 30 20 30 30	0001F. .
0 00C0:	30 20 30 30 20 36 38 20	30 30 20 30 30 20 30 30	0 00 68 00 00 00
0 00D0:	20 30 30 20 30 30 20 30	30 20 36 38 20 30 31 20	00 00 0 0 68 01
0 00E0:	30 30 20 30 30 20 30 30	20 30 30 20 30 30 20 2A	00 00 00 00 00 00
0 00F0:	20 30 78 30 30 30 30 32	30 2D 30 78 30 30 30 30	0x00002 0-0x000
0 0100:	32 46 0D 0A 30 30 20 30	30 20 30 30 20 30 30 20	2F..00 0 0 00 00
0 0110:	30 30 20 30 30 20 30 30	20 30 30 20 30 30 20 30	00 00 00 00 00 00
0 0120:	30 20 30 30 20 30 30 20	30 30 20 30 30 20 30 30	0 00 00 00 00 00 00
0 0130:	20 30 30 20 2A 20 30 78	30 30 30 30 33 30 2D 30	00 * 0x 000030-
0 0140:	78 30 30 30 30 33 46 0D	0A 30 31 20 30 33 20 30	x00003F. .01 03
0 0150:	30 20 30 30 20 30 30 20	30 30 20 30 30 20 30 30	0 00 00 00 00 00 00
0 0160:	20 30 30 20 66 30 20 65	30 20 30 30 20 30 66 20	00 f0 e 0 00 0f
0 0170:	66 30 20 30 30 20 30 30	20 2A 20 30 78 30 30 30	f0 00 00 * 0x00
0 0180:	30 34 30 2D 30 78 30 30	30 30 34 46 0D 0A 31 66	040-0x00 004F..1
0 0190:	20 30 35 20 30 30 20 30	30 20 30 30 20 30 30 20	05 00 0 0 00 00
0 01A0:	30 66 20 30 30 20 30 30	20 30 66 20 31 33 20 31	0f 00 00 0f 13
0 01B0:	61 20 33 63 20 33 34 20	30 31 20 30 30 20 2A 20	a 3c 34 01 00 *


STARTING POINT






PROGRAMMER AS VIRTUAL MACHINE



- Enabled retries
- Eliminate fear of bricking device
- Enabled exploration







PROGRAMMER AS VIRTUAL MACHINE

- > 12000 files
- Commercial software
- Proprietary software










Acrobat® Reader™ 5.0
There's more to Acrobat than Reader!

 **Kithara**

 **Chilkat Software**
inventing power with simplicity

 **BIOTRONIK**



RENESAS
Everywhere you imagine.




ATTACK SURFACE



ATTACK SURFACE



CHILKAT USED FOR ZIPPING AND ENCRYPTION



EXPLOIT FROM 2008 CONFIRMED OLD LIBRARY

Vulnerability Details : [CVE-2008-5002 \(1 Metasploit modules\)](#)

Insecure method vulnerability in the ChilkatCrypt2.ChilkatCrypt2.1 ActiveX control (ChilkatCrypt2.dll 4.3.2.1) in Chilkat Crypt ActiveX Component allows remote attackers to create and overwrite arbitrary files via the WriteFile method. NOTE: this could be leveraged for code execution by creating executable files in Startup folders or by accessing files using hcp:// URLs. NOTE: some of these details are obtained from third party information.

Publish Date : 2008-11-10 Last Update Date : 2017-09-28

[Collapse All](#) [Expand All](#) [Select](#) [Select&Copy](#)

[▼ Scroll To](#) [▼ Comments](#) [▼ External Links](#)

[Search Twitter](#) [Search YouTube](#) [Search Google](#)

– CVSS Scores & Vulnerability Types

CVSS Score 9.3

Confidentiality Impact Complete (There is total information disclosure, resulting in all system files being revealed.)

Integrity Impact Complete (There is a total compromise of system integrity. There is a complete loss of system protection, resulting in the entire system being compromised.)

Availability Impact Complete (There is a total shutdown of the affected resource. The attacker can render the resource completely unavailable.)

Access Complexity Medium (The access conditions are somewhat specialized. Some preconditions must be satisfied to exploit)

Authentication Not required (Authentication is not required to exploit the vulnerability.)

Gained Access None

Vulnerability Type(s) Execute Code

CWE ID [20](#)

CHILKAT SAYS NOT CHILKAT

```
Chilkat.Zip zip = new Chilkat.zip();
zip.UnlockComponent("License String");
zip.OpenZip = "filename.zip";
zip.get_Encryption(); // 0 = not encrypted
```

```
test@Linuxfjlas-MacBook-Pro ~/Downloads » python3 DetectEncryptionChilkat.py
This zip is NOT password-protected.
This zip is not encrypted.
Success!
Amount of files in zip: 317
test@Linuxfjlas-MacBook-Pro ~/Downloads »
```



Anders Been Wilhelmsen

to support ▾

Hello Chilkat support,

I have a perhaps unusual request. If possible, I would like to receive a (very) old version of the Chilkat Python library. I am not sure if it is possible to get an older python version if possible.

The reason for this is for our masters degree at NTNU we've stumbled over some files which apparently served as a check for the password. We would very much like to try and figure out what they were.

Any assistance would be greatly appreciated, thank you so much.

best regards,

Anders Been Wilhelmsen

MSc Communication Technology, NTNU

DEEP DIVE INTO CHILKAT

- How does Chilkat work?
 - Read documentation
 - Contacted Chilkat

Pages: [1]

Author

Topic: Old Chilkat version for legacy encrypted data (Read 3 times)

anderbw

Newbie



Posts: 0

Karma: +0/-0



Old Chilkat version for legacy encrypted data
« on: February 28, 2018, 07:46:48 AM »

Note: This message is awaiting approval by a moderator.

Hello there,

Come across some data files which we'd like to decrypt, they seem to be encrypted w
Would this be possible somehow? I'd prefer python on linux, but anything is acceptable

The reason I believe it uses the old version is that newer programs like zip or 7z does
saying "this zip has been encrypted with Chilkatsoft zip", plus a b64 string. Therefore

Also, could someone confirm what string is used for the b64 mode? If I understood rig
anyone could provide then, the following:

- A library to handle the encryption and decryption of this type of data
- The string that is encrypted as the password check
- What mode of operation was used in the legacy encryption modes, block size, type

Thank you so much for any help.

Pages: [1]



Hello anderbw

Show unread posts since last visit.

Show new replies to your posts.

March 20, 2018, 05:03:42 AM

Home Help Search Profile My Messages Members Logout

Chilkat Forum

An Error Has Occurred!

The topic or board you are looking for appears to be either missing or off limits to you.

Back

DEEP DIVE INTO CHILKAT

- How does Chilkat work?
 - Read documentation
 - Contacted Chilkat
 - Used Chilkat Forum

DEEP DIVE INTO CHILKAT

- How Chilkat should work

```
Chilkat.Zip zip = new Chilkat.zip();
zip.UnlockComponent("License String");
zip.Encryption = 3;
zip.EncryptionKeyLength = 256;
zip.EncryptPassword = "Password"
```

WHERE IS CHILKAT REFERENCED?


```
$ grep -Ri ChilkatZip2  
...  
Binary file Hcc.dll matches  
...
```

```
$ strings -a Execute/Hcc.dll| grep -i chilkat  
IChilkatZip2  
ChilkatZip2Lib_TLB\  
TChilkatZip2UnzipPercentDone  
TChilkatZip2WriteZipPercentDone  
TChilkatZip2FileUnzipped  
TChilkatZip2ToBeAdded  
TChilkatZip2FileAdded  
TChilkatZip2ToBeZipped  
TChilkatZip2FileZipped  
TChilkatZip2ToBeUnzipped  
TChilkatZip2  
TChilkatZip2  
ChilkatZip2Lib_TLB#  
ChilkatZip2Lib_TLB
```

DISASSEMBLING HCC.DLL

```
loc_463BC267:  
0A0 mov     edx, offset a9biotronikzip7 ; "9BiotronikZI  
0A0 mov     eax, [ebp+Zip20bj]  
0A0 call    ck_unlock  
0A0 mov     ecx, 3          ; encryption_mode_type  
0A0 mov     edx, 21h        ; encryption_property  
0A0 mov     eax, [ebp+Zip20bj] ; this  
0A0 call    set_property  
0A0 mov     ecx, 100h       ; keysize  
0A0 mov     edx, 22h        ; keysize_property  
0A0 mov     eax, [ebp+Zip20bj] ; this  
0A0 call    set_property  
0A0 mov     edx, offset aBiotronik ; "BIOTRONIK"  
0A0 mov     eax, [ebp+Zip20bj]  
0A0 call    setpw_or_encpw  
0A0 lea     eax, [ebp+var_4C]  
0A0 mov     edx, [ebp+var_8]  
0A0 call    @System@@WStrFromLStr$qqrr17System@WideStri  
0A0 mov     edx, [ebp+var_4C]  
0A0 mov     eax, [ebp+Zip20bj]  
0A0 call    sub_463BB960  
0A0 mov     eax, [ebp+var_4]  
0A0 cmp     dword ptr [eax+4], 0  
0A0 jnz     short loc_463BC2D6
```


DISASSEMBLING HCC.DLL



```
loc_463BC267:  
0A0 mov     edx, offset a9biotronikzip7 ; "9BiotronikZI  
0A0 mov     eax, [ebp+Zip20bj]  
0A0 call    ck_unlock  
0A0 mov     ecx, 3          ; encryption_mode_type  
0A0 mov     edx, 21h        ; encryption_property  
0A0 mov     eax, [ebp+Zip20bj] ; this  
0A0 call    set_property  
0A0 mov     ecx, 100h       ; keysize  
0A0 mov     edx, 22h        ; keysize_property  
0A0 mov     eax, [ebp+Zip20bj] ; this  
0A0 call    set_property  
0A0 mov     edx, offset aBiotronik ; "BIOTRONIK"  
0A0 mov     eax, [ebp+Zip20bj]  
0A0 call    setpw_or_encpw  
0A0 lea     eax, [ebp+var_4C]  
0A0 mov     edx, [ebp+var_8]  
0A0 call    @System@@WStrFromLStr$qqrr17System@WideStri  
0A0 mov     edx, [ebp+var_4C]  
0A0 mov     eax, [ebp+Zip20bj]  
0A0 call    sub_463BB960  
0A0 mov     eax, [ebp+var_4]  
0A0 cmp     dword ptr [eax+4], 0  
0A0 jnz     short loc_463BC2D6
```

```
Chilkat.Zip zip = new Chilkat.zip();  
zip.UnlockComponent("License String");  
zip.Encryption = 3;  
zip.EncryptionKeyLength = 256;  
zip.EncryptPassword = "Password"
```

- Have an attack vector
 - Can forge import files
- Can decrypt exported data
 - Including Marie's data
- Decrypted data exported in 2016
 - Still use Chilkat
 - Still use password-encrypted AES with «BIOTRONIK» as password



DEMO

S | Object Browser | ConsoleApp1 | Program.cs | C# | ConsoleApp1 | Main(string[] args) | Diagnostic Tools

Server Explorer | Toolbox

```
using System.Text;
using System.Threading.Tasks;
using Chilkat;

namespace ConsoleApp1
{
    class Program
    {
        static void Main(string[] args)
        {
            Chilkat.Zip zip = new Chilkat.Zip();
            zip.UnlockComponent("dummy_value");
            string filename = "C:\\Users\\IEUser\\Desktop\\2016_09_22_08_37_47_402_0.zip";
            zip.DecryptPassword = "BIOTRONIK";
            bool pw_ok = zip.VerifyPassword();
            if (pw_ok)
            {
                Console.WriteLine("Password is correct :-)");
            }
            else {
                Console.WriteLine(zip.LastErrorText);
            }

            zip.OpenZip(filename);
            zip.Extract("C:\\Users\\IEUser\\Desktop\\out_data");
        }
    }
}
```

100 % No issues found

Output
Show output from: Debug
'ConsoleApp1.exe' (CLR v4.0.30319: DefaultDomain): Loaded 'C:\Windows\Microsoft.Net\assembly\GAC_...
'ConsoleApp1.exe' (CLR v4.0.30319: DefaultDomain): Loaded 'C:\Users\IEUser\Source\Repos\ConsoleAp...
'ConsoleApp1.exe' (CLR v4.0.30319: ConsoleApp1.exe): Loaded 'C:\Users\IEUser\Source\Repos\Console...
The program '[5808] ConsoleApp1.exe' has exited with code 0 (0x0).

Diagnostic Tools
Diagnostics session: 4 seconds (4.591 s selected)
Events
Process Memory (MB)
10 10

Open
Extract All...
Open with...
Share with
Restore previous versions
Send to
Cut
Copy
Create shortcut
Delete
Rename
Properties

Open in new window

Extract All... (highlighted)

Open with...

Share with

Restore previous versions

Send to

Cut

Copy

Create shortcut

Delete

Rename

Properties

2016_09_22_08_37_47_402_0.zip Compressed (zipped) Folder

vs_community_1738190835.156218 3892.exe Visual Studio Installer

2016_09_22_08_37_47_402_0.zip Date modified: 7/15/2019 5:12 AM
Compressed (zipped) Folder Size: 558 KB

Date created: 7/15/2019 5:14 AM



DISCLOSURE PROCESS

- We have informed Biotronik of our findings.
- Biotronik has confirmed that they will update the import/export function
- Biotronik claims none of these vulnerabilities represents an uncontrolled risk as defined by FDA

- Thesis download link: <https://anderbw.github.io>

QUESTIONS?

Thesis download link: <https://anderbw.github.io>



THANK YOU!

ANDERS B.WILHELMSEN @anderbw

EIVIND S. KRISTIANSEN @Skjelmo

MARIE MOE @MarieGMoe

Thesis download link: <https://anderbw.github.io>